



Scabil

Interchangeable handpiece allows multiple beam patterns for various applications and treatments.



Fractional & Surgical mode

- Equipped the latest RF metal tube provides a uniform and perfect beam pattern.
Excellent treatment effect, Short Down time, Reduce pain
- Dual Hand-piece (Pencil & Fractional type) and Dual Mode (Fractional & Surgical mode) load

Indication

- Acne Scars
- Skin resurfacing
- Pore minimization
- Reduce Wrinkles
- Skin tone & Texture
- Pigmented Skin
- Stretch mark
- Surgical scars
- Burn Injuries

Scanning patterns



Specifications

Irradiation Type	Diode laser
Wave length	808nm
Repetition Rate	Single, up to 1 ~ 20Hz
Pulse Duration	10 -800ms
Maximum Energy	1J/cm ² ~ 233J/cm ²

WON TECH Co., Ltd.

www.wtlaser.com

Address : Pangyo-ro 242, Bundang-gu, Seongnam-si, Gyeonggi-do, Korea
(Sampyeong-dong, A-601 Pangyo digital contents park)

Tel : +82.70.7882.8657 **Fax** : +82.70.7552.8658

Email : asia@wtlaser.com / middleeast@wtlaser.com / europa@wtlaser.com / china@wtlaser.com / america@wtlaser.com

www.wtlaser.com
Scabil

