

# 6D TAPE<sup>®</sup>

Six-Dimensional Method

## Premium KINESIOLOGY TAPE WITH HANDLES



**Its main feature is to mobilize tissue and muscles.**

**Is designed for physio and lymph therapy, as well as self-administered care.**

The therapist or the client can activate and mobilize the tissue in different directions with 6D Tape. The mobilization has a positive impact on the lymphatic system, blood circulation, fascia, muscles and joints without limiting their functions.

There are plenty of variations of 6D taping techniques depending on desirable results and the object of taping.

Feel free to use a combination of traditional kinesiology taping techniques with 6D Tape.

The product is an effective and affordable tool for professionals. With minimal training, customers / patients can continue with the active care at home in order to increase the speed of recovery. The 6D Tape remains in place on the skin from 3 to 10 days depending on the intensity and quantity of mobilization

### Unique 6D Handles

The unique handles are designed to mobilize the tissue in 6 dimensions.

### Material

This base material is unique and skin friendly kinesiology tape. The kinesiology tape is made of 98% cotton. The material is breathable and water-resistant. The base material is elastic, with 135%-170% stretch ability.

The handle material is non-elastic.

### Adhesive

Soft gel adhesive reduces pain and damage to skin on removal by conforming to the rough surface of skin.

The acrylic adhesive is heat-activated with a gentle rub, and the bond reaches its maximum strength in 15 minutes time. Moreover, it realizes the high adhesive with the conformity to the skin's rough surface.

FROM PASSIVE TAPING

TO ACTIVE TREATMENT



PATENT PENDING

Exclusive Distributor :

**ABSOLUTE MS (S) PTE LTD**

73 Ubi Road 1 #07-48 Oxley Bizhub Singapore 408733

Tel: 6509 0091

[www.absolutems.com.sg](http://www.absolutems.com.sg)